



Versatility When you Need It Most.

Terumo Medical Corporation is the preferred NuEdge CONNECT program for Hypodermic products.

Did you know that pairing Terumo's **Venoject** Multi-Sample Luer Adapter with Terumo's **Surshield™** Safety Winged Infusion sets transform a standard infusion set into a winged blood collection set? Simply connect Terumo's **Venoject** Multi-Sample Luer Adapter to any of Terumo's **Surshield™** Safety Winged Infusion sets to create a safety winged blood collection set. The **Venoject** Multi-sample Luer adapter is sterile and disposable, allowing multiple draws into evacuated tubes and accept standard tube holders.



Terumo's Surshield™ Safety Winged Infusion Sets are unique from other safety winged butterflies in that most safety winged butterflies are designed to be activated one way and require a second hand; Surshield™ provides **two** different one-handed activation options:

- Hard Surface
- Thumb



**Refer to product labels and package insert for complete warnings, precautions, potential complications, and instructions for use*

Two one-handed activation techniques give options without sacrificing personal technique. Terumo designs injection and infusion therapy devices that help clinicians consistently perform at their very best.

When versatility and superior quality is of the utmost importance, look no further than Terumo's Surshield™ and Luer Adapter combination for your practice.

To view **Terumo's Surshield™** Safety Winged Infusion Sets IFU video please visit:

<https://www.terumotmp.com/products/infusion-and-blood-collection/surshield-safety-winged-infusion-sets.html>

PM-05422