

# HAEVOLANCE<sup>®</sup> Plus

Safety, efficacy and comfort in two simple steps

CE 0344

## Features and benefits

### SAFETY

- Single-use only
- automatic needle retraction and lock out
- needlestick injury and cross-infection prevention

### EFFICACY

- wide product range – available in 6 different sizes
- derives proper amount of blood for different capillary tests
- dedicated versions for heel stick blood tests on neonates

### COMFORT

- easy to use in just two steps
- design to control needle depth



# HAEMOLANCE<sup>+</sup> Plus

Safety, efficacy and comfort in two simple steps



## Paediatric

1,5 mm blade  
1,2 mm depth

Cat. no. 7375 / 100ct box  
Cat. no. 7365 / 200ct box



## Max Flow

1,5 mm blade  
1,6 mm depth

Cat. no. 7374 / 100ct box  
Cat. no. 7364 / 200ct box



## High Flow

18G needle  
1,8 mm depth

Cat. no. 7373 / 100ct box  
Cat. no. 7363 / 200ct box



## Normal Flow

21G needle  
1,8 mm depth

Cat. no. 7372 / 100ct box  
Cat. no. 7362 / 200ct box



## Low Flow

25G needle  
1,4 mm depth

Cat. no. 7371 / 100ct box  
Cat. no. 7361 / 200ct box



## Micro Flow

28G needle  
1,6 mm depth

Cat. no. 7370 / 100ct box  
Cat. no. 7360 / 200ct box

## INSTRUCTIONS FOR USE:



Twist off the protective cap, and then pull it straight out.



Place HAEMOLANCE PLUS<sup>®</sup> on the desired puncture site and press the button to activate the device.



Gently apply intermittent pressure near the puncture site to obtain the required blood volume.